

## MODULAR FLASH POINT TESTER

**MULTIFLASH AUTOMATIC MULTI METHOD FLASH POINT TESTER  
(PENSKY-MARTENS, ABEL, TAG, CLEVELAND & SMALL SCALE)**

Multiflash is a universal fully automated instrument designed for testing the flash point of fuels, lubricating oils, solvents, chemicals, waste materials and many other substances.

The instrument comprises a universal base unit that houses the display, control and data acquisition functions, and a series of interchangeable modules that allow the implementation of the required standard test methods.

Multiflash is ruggedly designed to provide simple and reliable operation and is specifically engineered to ensure a very low cost of maintenance.



BIO-FUEL TESTING



CALIBRATION &amp; VERIFICATION MATERIALS



SETALERT COMPATIBLE



TECHNICAL DATASHEETS AVAILABLE

**PENSKY-MARTENS (A, B & C Closed Cup)**

ASTM D93; E502; IP 34; IP 404 (obs);  
BS 2000 Parts 34 & 404 (obs); ISO 15267;  
EN ISO 2719; JIS K2265; NF M07-019;  
DIN 51 758; FTM 791 1102; EPA 1010A;  
GB/T 261

**ABEL (Closed Cup and Equilibrium)**


IP 170; IP 304 (obs) Part 2; IP 492; BS 2000  
Part 170; BS 3900 Parts A8 and A9; BS  
6664 Parts 1 and 2; EN ISO 13736; ISO 1523

**CLEVELAND (Open Cup)**

ASTM D92; IP 36; EN ISO 2592;  
BS 2000 Parts 36 & 403 (obs); DIN 51 376;  
NF T60-119; JIS K2265; AASHTO T48

**TAG (Closed Cup)**

ASTM D56; D3934; D3941; E502;  
IP 304 (obs) Parts 1 & 2; IP 491; IP 492;  
ISO 1516; ISO 1523; FTM 791 1101;  
JIS K2265

**SMALL SCALE Rapid Equilibrium  
(Setaflash) Closed Cup**

ASTM D3278; D3828; D7236; E502;  
IP 303 (obs); IP 523; IP 524;  
IP 534; ISO 3679; ISO 3680; BS EN 456;  
BS 6664 Parts 3 & 4; BS 3900 Part A13;  
UN Class 3 Non-viscous Flammable  
Liquids; CHIPS Regulations;  
Classification of Dangerous Goods for  
Carriage: Viscous and Non-viscous  
Liquids; EPA 1020 A & B



34000-0 UNIVERSAL BASE UNIT

## MODULAR FLASH POINT TESTER

- **ALL MAJOR TEST METHODS**

Pensky-Martens (A,B & C Closed Cup and Chinese Method: GB)  
Abel (Closed Cup & Equilibrium)  
TAG (Closed Cup & Equilibrium)  
Cleveland (Open Cup)  
Small Scale Rapid Equilibrium (Setaflash)

- **UNIVERSAL OPERATION**

Unique '5 in 1' modular solution for flash point testing  
All primary flash point tests from one instrument  
Universal base, interchangeable test modules  
Gas or electric ignition

- **FULLY AUTOMATIC**

Automatic set up of test parameters  
Automatic heating control, ignition & flash detection  
Automatic test results & end of test warning  
Automatic barometric correction

- **UNRIVALLED EASE OF USE**

Minimum operator input  
Integral self diagnosis

- **LOW COST**

Modular design optimises capital outlay  
Instrument can be configured to match the requirements of the user  
Simple to extend range of test methods by adding modules

- **STRICT CONFORMITY TO TEST METHODS**

Precise results  
Procedures in conformance with method of test  
Critical components designed in exact conformance with test methods  
Easy and quick calibration and verification procedures

- **FAST & RELIABLE PERFORMANCE**

Designed for long service life  
Reliable and proven technology  
Simple operation  
Error messaging  
Calibration & diagnostics mode

- **IN-BUILT SAFETY**

3 independent hardware safety systems  
Fire sensors, over-temperature protection and automatic gas shut off  
Automatic sample, bath fluid and probe checks prior to heating  
Integral software limits  
Suitable for unattended operation  
Pensky-Martens D93 10°C safety protocol

- **RANGE OF ACCESSORIES**

Clear draught screens  
Printer  
Remote alarm system  
Rapid cooling modules & cryostats  
Calibration and verification instruments and materials



## MODULAR FLASH POINT TESTER

## 34100-2 MULTIFLASH - PENSKY-MARTENS CLOSED CUP MODULE

ASTM D93; E502; IP 34; IP 404 (obs); BS 2000 Parts 34 & 404 (obs); ISO 15267; EN ISO 2719; JIS K2265; NF M07-019; DIN 51 758; FTM 791 1102; EPA 1010A; GB/T 261

- Procedures A, B and C
- Flexible automated search mode
- Ambient to 400°C temperature range
- Gas or electric ignition
- Independent fire detector

The Pensky-Martens Closed Cup Module comprises a Test Area, DIPS pod (Dipping, Ignition, Pilot, Stirrer), flash detector, sample temperature probe, and all necessary cables to interconnect with the 34000-0 Universal Base Unit.

The Test Area consists of an electrically heated solid block/air bath with a removeable sample cup and lid. The lid is equipped with a stirrer, rotating shutter, and ports for the flash detector and sample temperature probe.

The DIPS pod contains a fire detection probe, stirrer drive, shutter actuator and ignitor dipping mechanism with a coaxial gas jet/electric hot wire ignitor.

The sample temperature probe is a class A PRT, and flash detection is by a thermocouple.

Safety features include an initial dip to check that the sample is below its flash point, a PRT check to ensure that the correct probe is connected and is immersed in the sample, and new 10°C dipping protocol. When in search mode, the software limits the maximum temperature to 50°C above "expected flash". A bath temperature probe protects against over-temperature, and an independent fire detection probe is mounted over the test area.

The optional 34008-0 Cooling Boost Module allows the test area to be rapidly cooled in between tests by a blast of compressed air.

**SUPPLIED WITH:** DIPS pod, Test Area, Ignitor, Flash Detector, Sample PRT, interconnection leads, and instruction manual.

*Note: 34000-0 Universal Base Unit required for operation.*



34100-2 mounted on a 34000-0 Base Unit with a 34010-0 Draught Screen fitted



BIO-FUEL TESTING



CALIBRATION & VERIFICATION MATERIALS



SETALERT COMPATIBLE



TECHNICAL DATASHEETS AVAILABLE

| SPECIFICATIONS                                  |  |
|---|--|
| <b>Method:</b>                                  | Pensky-Martens Closed Cup  |
| <b>Tests:</b>                                   | P-M Procedures A, B & C<br>GB/T 261<br>Search mode for unknown samples                                   |
| <b>Temperature Range:</b>                       | Ambient to 400°C   |
| <b>Barometric Correction:</b>                   | Automatic  |
| <b>Bath:</b>                                    | Solid Block/Air  |
| <b>Heating:</b>                                 | Bottom and side heaters  |
| <b>Post Test Cooling:</b>                       | Fan assisted or<br>Forced Air (Requires 34008-0 Cooling Boost Module option)                             |
| <b>Post Test Cool-down Time (90°C to 50°C):</b> | 36 minutes (Fan)<br>6 minutes (Forced Air)   |
| <b>Ignitor:</b>                                 | Quick-fit coaxial gas jet/electric hot wire, user selectable (Electric hot wire pilot light for gas jet) |
| <b>Temperature Probes:</b>                      | Class A Platinum Resistance Thermometer, °C or °F  |
| <b>Flash Detection:</b>                         | Thermocouple   |
| <b>Stirrer Speed:</b>                           | 105 & 250 rpm  |
| <b>Gas Supply (gas jet ignitor):</b>            | Laboratory gas, 3kPa (0.44psi) maximum pressure  |
| <b>Air Supply (forced air cooling):</b>         | Dry compressed air, 670kPa (100psi) maximum pressure   |
| <b>Voltage:</b>                                 | 110/120V or 220/240V, 50/60Hz, (switchable)  |
| <b>Power:</b>                                   | 1kW  |
| <b>Size (HxWxD):</b>                            | 40 x 45 x 34cm   |
| <b>Weight:</b>                                  | 20kg   |

## ACCESSORIES - Pensky-Martens Closed Cup Module

## OPTIONAL

|         |   |
|---------|---|
| 00100-0 | SETALERT REMOTE PAGING SYSTEM                 |
| 23680-0 | SPIRIT LEVEL                                  |
| 34003-0 | RS232 CABLE, for connecting instrument to PC. |
| 34008-0 | COOLING BOOST MODULE                          |
| 34010-0 | DRAUGHT SCREEN (recommended for gas tests).   |
| 81002-2 | PRINTER                                       |

## CALIBRATION VERIFICATION &amp; TRAINING

|           |  |
|-----------|--|
| 60005-002 | DRY WELL PROBE CALIBRATOR, (for details see page 144).                                     |
| 99851-0   | MULTI-TEST VERIFICATION MATERIAL, D93 - IP 34 Gas Oil, (for details see page 150).         |
| 99853-0   | MULTI-TEST VERIFICATION MATERIAL, D93 - IP 34 Lubricating Oil, (for details see page 151). |
| 99908-0   | BIO MULTI-TEST VERIFICATION MATERIAL, Gas Oil, (for details see page 151).                 |
| 99880-0   | FLASHCHECK, secondary verification material, (for details see page 156).                   |
| 99910-0   | PORTABLE DIGITAL BAROMETER, 750 to 1500 mbar absolute, (for details see page 145).         |
| 99960-0   | DIGITAL PHOTO TACHOMETER, (for details see page 147).                                      |

## CONSUMABLES - Pensky-Martens Closed Cup Module

| 2 Years/2000 Tests |                                 | Qty |
|--------------------|---------------------------------|-----|
| 34100-001          | FLASH DETECTOR                  | 2   |
| 34100-002          | SAMPLE PRT PROBE                | 1   |
| 34100-004          | CO-AXIAL IGNITOR                | 2   |
| 81002-001          | PRINTER PAPER, pack of 20 rolls | 1   |

**34200-0 MULTIFLASH - ABEL CLOSED CUP MODULE**

IP 170; IP 304 (obs) Part 2; IP 492; BS 2000 Part 170; BS 3900 Parts A8 and A9; BS 6664 Parts 1 & 2; EN ISO 13736; ISO 1523

- Equilibrium tests
- Non-Equilibrium tests
- Flash/No-Flash tests
- Flexible automated search mode
- Ambient to 93°C (5°C to 93°C with optional cryostat)
- Gas or electric ignition

The Abel Closed Cup Module comprises a Test Area, DIPS pod (Dipping, Ignition, Pilot, Stirrer), flash detector, sample temperature probe, outer and inner bath temperature probes and all necessary cables to interconnect with the 34000-0 Universal Base Unit.

The Test Area consists of an outer liquid filled bath surrounding an inner liquid filled bath and a removeable sample cup and lid. The outer bath is heated by the hotplate within the 34000-0 Universal Base Unit, has a port for a temperature probe and connections for a 34008-0 Cooling Boost Module, allowing circulation of cooling fluid via a cryostat for sub-ambient tests. The inner bath also has provision for a temperature probe. The lid of the sample cup is equipped with a stirrer, sliding shutter, and ports for the flash detector and sample temperature probe.

The DIPS pod contains a stirrer drive, shutter actuator and ignitor dipping mechanism with a coaxial gas jet/electric hot wire ignitor.

The outer bath, inner bath and sample temperature probes are class A PRTs, and flash detection is by a thermocouple.

Safety features include an initial dip to check that the sample is below its flash point at the start of the test and PRT checks to ensure that the correct probes are connected and are immersed in bath liquids or sample. When in search mode, the software limits the maximum temperature to 50°C above "expected flash". An over-temperature cut-out is mounted on the base flange of the outer bath.

**SUPPLIED WITH:** DIPS pod, Test Area, Ignitor, Flash Detector, Inner and Outer Bath PRTs, Sample PRT, interconnection leads, instruction manual.

*Note: 34000-0 Universal Base Unit and 34008-0 Cooling Boost Module are required for operation. Sub-ambient tests down to 5°C require a 34006-0 Cryostat.*

| SPECIFICATIONS                                  |  |
|---|--|
| <b>Method:</b>                                  | Abel Closed Cup  |
| <b>Tests:</b>                                   | Equilibrium<br>Non-Equilibrium<br>Flash/No Flash   |
| <b>Temperature Range:</b>                       | Ambient to 93°C (5°C to 93°C with optional cryostat)   |
| <b>Barometric Correction:</b>                   | Automatic  |
| <b>Bath:</b>                                    | Liquid/Liquid  |
| <b>Heating:</b>                                 | Hot Plate  |
| <b>Post Test Cooling:</b>                       | Fan assisted or<br>Circulating liquid/external cryostat (optional)                                       |
| <b>Post Test Cool-down Time (90°C to 50°C):</b> | 45 minutes (Fan)<br>3 minutes (external cryostat)  |
| <b>Ignitor:</b>                                 | Quick-fit coaxial gas jet/electric hot wire, user selectable (Electric hot wire pilot light for gas jet) |
| <b>Temperature Probes:</b>                      | Class A Platinum Resistance Thermometers   |
| <b>Flash Detection:</b>                         | Thermocouple   |
| <b>Stirrer Speed:</b>                           | 30 rpm   |
| <b>Gas Supply (gas jet ignitor):</b>            | Laboratory gas, 3kPa (0.44psi) maximum pressure  |
| <b>Cooling Supply:</b>                          | Water (5°C and above)<br>Water/Ethylene Glycol 50/50 (below 5°C)<br>125kPa (18psi) maximum pressure      |
| <b>Voltage:</b>                                 | 110/120V or 220/240V, 50/60Hz, (switchable)  |
| <b>Power:</b>                                   | 600W   |
| <b>Size (HxWxD):</b>                            | 47 x 45 x 34cm   |
| <b>Weight:</b>                                  | 13kg   |



34200-0 mounted on a 34000-0 Base Unit with a 34012-0 Draught Screen fitted

-  CALIBRATION & VERIFICATION MATERIALS
-  SETALERT COMPATIBLE
-  TECHNICAL DATASHEETS AVAILABLE

**ACCESSORIES - Abel Closed Cup Module****REQUIRED**

34008-0 COOLING BOOST MODULE

**OPTIONAL**

00100-0 SETALERT REMOTE PAGING SYSTEM  
23680-0 SPIRIT LEVEL  
34003-0 RS232 CABLE, for connecting instrument to PC.  
34006-0 CRYOSTAT (-20°C).  
34007-0 SHUTTER LUBE  
34012-0 DRAUGHT SCREEN (recommended for gas tests).  
81002-2 PRINTER

**CALIBRATION VERIFICATION & TRAINING**

60005-002 DRY WELL PROBE CALIBRATOR, (for details see page 144).  
99850-0 MULTI-TEST VERIFICATION MATERIAL, IP 170 Kerosine, (for details see page 150).  
99910-0 PORTABLE DIGITAL BAROMETER, 750 to 1500 mbar absolute, (for details see page 145).  
99960-0 DIGITAL PHOTO TACHOMETER, (for details see page 147).

**CONSUMABLES - Abel Closed Cup Module**

2 Years/2000 Tests

|   | Qty |
|---|-----|
| 34100-002 SAMPLE PRT PROBE                | 1   |
| 34200-003 OUTER BATH PRT PROBE            | 1   |
| 34200-005 CO-AXIAL IGNITOR                | 2   |
| 34700-001 FLASH DETECTOR                  | 2   |
| 81002-001 PRINTER PAPER, pack of 20 rolls | 1   |

## MODULAR FLASH POINT TESTER

## 34300-0 MULTIFLASH - CLEVELAND OPEN CUP MODULE

ASTM D92; IP 36; EN ISO 2592; BS 2000 Parts 36 & 403 (obs);  
DIN 51 376; NF T60-119; JIS K2265; AASHTO T48

- Flash Point & Fire Point
- Flexible automated search mode
- Ambient to 400°C temperature range
- Gas or electric ignition
- Automatic snuffer
- Independent fire detector

The Cleveland Open Cup Module comprises a Test Area, DIPS pod (Dipping, Ignition, Pilot, Stirrer), flash/fire detector, sample temperature probe, and all necessary cables to interconnect with the 34000-0 Universal Base Unit.

The Test Area consists of a removeable sample cup in a thermal insulation ring, and a probe block for mounting the flash/fire detector and sample temperature probe. The sample cup is heated by the hotplate within the 34000-0 Universal Base Unit.

The DIPS pod contains the ignitor sweeping mechanism with a coaxial gas jet/electric hot wire ignitor, and a snuffer that extinguishes any flames in the sample cup at the end of the test.


The sample temperature probe is a class A PRT, and flash/fire detection is by an ionisation ring.

Safety features include an initial sweep to check that the sample is below its flash point at the start of the test and a PRT check to ensure that the correct probe is connected and is immersed in sample. When in search mode, the software limits the maximum temperature to 50°C above "expected flash". The snuffer automatically deploys at the end of the test or if the independent fire detection probe (mounted over the test area) detects a flame. The hotplate in the 34000-0 is equipped with a cut-out to protect against over-temperature.

The optional 34008-0 Cooling Boost Module allows the test area to be rapidly cooled in between tests by a blast of compressed air.

**SUPPLIED WITH:** DIPS pod, Test Area, Ignitor, Flash/Fire Detector, Sample PRT, interconnection leads, and instruction manual.

*Note: 34000-0 Universal Base Unit required for operation.*



34300-0 mounted on a 34000-0 Base Unit with a 34012-0 Draught Screen fitted



CALIBRATION & VERIFICATION MATERIALS



SETALERT COMPATIBLE



TECHNICAL DATASHEETS AVAILABLE

| SPECIFICATIONS                              |   |
|---|---|
| Method:                                     | Cleveland Open Cup  |
| Tests:                                      | Flash Point and Fire Point<br>Non-Equilibrium<br>Search mode for unknown samples                            |
| Temperature Range:                          | Ambient to 400°C  |
| Barometric Correction:                      | Automatic   |
| Bath:                                       | N/A   |
| Heating:                                    | Hotplate  |
| Post Test Cooling:                          | Fan assisted or<br>Forced Air (Requires 34008-0 Cooling Boost Module option)                                |
| Post Test Cool-down Time<br>(90°C to 50°C): | 9 minutes (Fan)<br>6 minutes (Forced Air)   |
| Ignitor:                                    | Quick-fit coaxial gas jet/electric hot wire, user selectable<br>(Electric hot wire pilot light for gas jet) |
| Temperature Probe:                          | Class A Platinum Resistance Thermometer, °C or °F   |
| Flash/Fire Detection:                       | Ionisation Ring   |
| Gas Supply (gas jet ignitor):               | Laboratory gas, 3kPa (0.44psi) maximum pressure   |
| Air Supply (forced air cooling):            | Dry compressed air, 670kPa (100psi) maximum pressure  |
| Voltage:                                    | 110/120V or 220/240V, 50/60Hz, (switchable)   |
| Power:                                      | 1kW   |
| Size (HxWxD):                               | 34 x 45 x 42cm  |
| Weight:                                     | 20kg  |

## ACCESSORIES - Cleveland Open Cup Module

## OPTIONAL

|         |   |
|---------|---|
| 00100-0 | SETALERT REMOTE PAGING SYSTEM                 |
| 23680-0 | SPIRIT LEVEL                                  |
| 34003-0 | RS232 CABLE, for connecting instrument to PC. |
| 34008-0 | COOLING BOOST MODULE                          |
| 34012-0 | DRAUGHT SCREEN (recommended for gas tests).   |
| 34301-0 | HUSSAMI KIT                                   |
| 81002-2 | PRINTER                                       |

## CALIBRATION VERIFICATION &amp; TRAINING

|           |  |
|-----------|--|
| 60005-002 | DRY WELL PROBE CALIBRATOR, (for details see page 144).                             |
| 99880-0   | FLASHCHECK, secondary verification material, (for details see page 156).           |
| 99910-0   | PORTABLE DIGITAL BAROMETER, 750 to 1500 mbar absolute, (for details see page 145). |

## CONSUMABLES - Cleveland Open Cup Module

| 2 Years/2000 Tests |                                       | Qty |
|--------------------|---------------------------------------|-----|
| 34100-004          | CO-AXIAL IGNITOR .....                | .2  |
| 34300-001          | IONISATION RING .....                 | .2  |
| 34300-002          | SAMPLE PRT PROBE .....                | .1  |
| 81002-001          | PRINTER PAPER, pack of 20 rolls ..... | .1  |

**34400-0 TAG CLOSED CUP MODULE**

ASTM D56; D3934; D3941; E502; IP 304 (obs) Part 2;  
IP 491; IP 492; ISO 1516; ISO 1523; FTM 791 1101; JIS K2265

- Equilibrium tests
- Non-Equilibrium tests
- Flash/No-Flash tests
- Flexible automated search mode
- Ambient to 93°C (5°C to 93°C with optional cryostat)
- Gas or electric ignition

The TAG Closed Cup Module comprises a Test Area, DIPS pod (Dipping, Ignition, Pilot, Stirrer), flash detector, sample temperature probe, bath temperature probe and all necessary cables and pipes to interconnect with the 34000-0 Universal Base Unit and 34008 Cooling Boost Module.

The Test Area consists of a liquid filled bath with an internal cooling coil, a removeable sample cup and lid, and an overflow vessel. The bath is heated by the hotplate within the 34000-0 Universal Base Unit. The cooling coil is connected to the 34008-0 Cooling Boost Module by flexible pipes allowing circulation of water or cooling fluid via a cryostat for sub-ambient tests. The top of the bath has a port for a temperature probe. The lid of the sample cup is equipped with a sliding shutter, and ports for the flash detector and sample temperature probe.

The DIPS pod contains the shutter actuator and ignitor dipping mechanism with a coaxial gas jet/electric hot wire ignitor.

The sample and bath temperature probes are class A PRTs, and flash detection is by a thermocouple.

Safety features include an initial dip to check that the sample is below its flash point at the start of the test and PRT checks to ensure that the correct probes are connected and are immersed in bath liquid or sample. When in search mode, the software limits the maximum temperature to 50°C above "expected flash". An over-temperature cut-out is mounted on the base flange of the bath.

**SUPPLIED WITH:** DIPS pod, Test Area, Ignitor, Flash Detector, Bath PRT, Sample PRT, interconnection leads and instruction manual.

*Note: Requires 34000-0 Universal Base Unit and 34008-0 Cooling Boost Module. For sub-ambient tests down to 5°C a Cryostat (34006-0) is also required.*

| SPECIFICATIONS                                      |   |
|---|---|
| <b>Method:</b>                                      | TAG Closed Cup  |
| <b>Tests:</b>                                       | Non-Equilibrium<br>Equilibrium<br>Flash/No Flash  |
| <b>Temperature Range:</b>                           | Ambient to 93°C<br>(5°C to 93°C with optional cryostat)   |
| <b>Barometric Correction:</b>                       | Automatic   |
| <b>Bath:</b>  | Liquid  |
| <b>Heating:</b>                                     | Hot Plate   |
| <b>Post Test Cooling:</b>                           | Fan assisted or<br>Circulating liquid/external cryostat (Requires 34008-0<br>Cooling Boost Module option)   |
| <b>Post Test Cool-down Time<br/>(90°C to 50°C):</b> | 40 minutes (Fan)<br>3 minutes (external cryostat)   |
| <b>Ignitor:</b>                                     | Quick-fit coaxial gas jet/electric hot wire, user selectable<br>(Electric hot wire pilot light for gas jet) |
| <b>Temperature Probes:</b>                          | Class A Platinum Resistance Thermometers  |
| <b>Flash Detection:</b>                             | Thermocouple  |
| <b>Gas Supply (gas jet ignitor):</b>                | Laboratory gas, 3kPa (0.44psi) maximum pressure   |
| <b>Cooling Supply:</b>                              | Water (5°C and above)<br>Water/Ethylene Glycol 50/50 (below 5°C)<br>125kPa (18psi) maximum pressure         |
| <b>Voltage:</b>                                     | 110/120V or 220/240V, 50/60Hz, (switchable)   |
| <b>Power:</b>                                       | 1kW   |
| <b>Size (HxWxD):</b>                                | 42 x 45 x 34cm  |
| <b>Weight:</b>                                      | 20kg  |



34400-0 mounted on 34000-0 with 34011-0 Draught Screen



CALIBRATION & VERIFICATION MATERIALS



SETALERT COMPATIBLE



TECHNICAL DATASHEETS AVAILABLE

**ACCESSORIES - TAG Closed Cup Module****REQUIRED**

34008-0 COOLING BOOST MODULE

**OPTIONAL**

00100-0 SETALERT REMOTE PAGING SYSTEM  
23680-0 SPIRIT LEVEL  
34003-0 RS232 CABLE, for connecting instrument to PC.  
34006-0 CRYOSTAT (-20°C).  
34007-0 SHUTTER LUBE  
34011-0 DRAUGHT SCREEN (recommended for gas tests).  
81002-2 PRINTER

**CALIBRATION VERIFICATION & TRAINING**

60005-002 DRY WELL PROBE CALIBRATOR, (for details see page 147).  
99880-0 FLASHCHECK, secondary verification material, (for details see page 156).  
99910-0 PORTABLE DIGITAL BAROMETER, 750 to 1500 mbar absolute, (for details see page 145).

**CONSUMABLES - TAG Closed Cup Module****2 Years/2000 Tests**

|   | Qty |
|---|-----|
| 34100-001 FLASH DETECTOR .....                  | 2   |
| 34400-002 SAMPLE PRT PROBE .....                | 1   |
| 34400-003 BATH PRT PROBE .....                  | 1   |
| 34400-005 CO-AXIAL IGNITOR .....                | 2   |
| 81002-001 PRINTER PAPER, pack of 20 rolls ..... | 1   |

## MODULAR FLASH POINT TESTER

## 34700-0 SMALL SCALE (SETAFLASH) CLOSED CUP MODULE

ASTM D3278; D3828; D7236; E502; IP 303 (obs) Parts 1 & 2; IP 523; IP 524; IP 534; ISO 3679; ISO 3680; BS 6664 Parts 3 & 4; BS 3900 Part A13; BS EN 456; EPA 1020 A & B; UN Class 3 Non-viscous Flammable Liquids; CHIPS Regulations; Classification of Dangerous Goods for Carriage: Viscous and Non-viscous Liquids

- Rapid Equilibrium Tests
- Non-Equilibrium tests
- Flash/No-Flash tests
- Flexible automated search (ramp) mode
- Ambient to 300°C (-10°C to 300°C with optional cryostat)
- Gas or electric ignition

The Small Scale (Setaflash) Closed Cup Module comprises a Test Area, DIPS pod (Dipping, Ignition, Pilot, Stirrer), flash detector, and all necessary cables to interconnect with the 34000-0 Universal Base Unit.

The Test Area consists of an electrically heated solid block bath/sample cup with an integral temperature probe. The bath has a network of internal passages and can be connected to a 34008-0 Cooling Boost Module, allowing circulation of cooling fluid via a cryostat for sub-ambient tests. The hinged lid is equipped with a shutter and ports for the flash detector and injecting sample.

The DIPS pod contains the shutter actuator and ignitor dipping mechanism with a coaxial gas jet/electric hot wire ignitor.

The temperature probe is a class A PRT and flash detection is by a precision thermocouple.

Safety features include a check to ensure that the sample probe and flash detector are connected. When in search mode, the software limits the maximum temperature to 50°C above "expected flash".

**SUPPLIED WITH:** DIPS pod, Test Area, Ignitor, Flash Detector, Sample PRT, interconnection leads, and instruction manual.

*Note: 34000-0 Universal Base Unit required for operation. For sub-ambient tests a Cooling Boost Module (34008-0) and Cryostat (34006-0) are also required.*



34700-0 mounted on a 34000-0 Base Unit with a 34010-0 Draught Screen fitted



BIO-FUEL TESTING



CALIBRATION & VERIFICATION MATERIALS



SETALERT COMPATIBLE



TECHNICAL DATASHEETS AVAILABLE

## SPECIFICATIONS

|   |   |
|---|---|
| <b>Method:</b>                                  | Small Scale Closed Cup  |
| <b>Tests:</b>                                   | Rapid Equilibrium<br>Non-Equilibrium<br>Flash/No Flash<br>Search mode for unknown samples   |
| <b>Temperature Range:</b>                       | Ambient to 300°C (-10°C to 300°C with optional cryostat)  |
| <b>Barometric Correction:</b>                   | Automatic   |
| <b>Bath:</b>                                    | Solid Block   |
| <b>Heating:</b>                                 | Side heaters  |
| <b>Post Test Cooling:</b>                       | Fan assisted<br>Forced air or circulating liquid option (Requires 34008-0 Cooling Boost Module and 34006-0 Cryostat)  |
| <b>Post Test Cool-down Time (90°C to 50°C):</b> | 36 minutes (Fan)<br>2 minutes (Cryostat)  |
| <b>Ignitor:</b>                                 | Quick-fit coaxial gas jet/electric hot wire, user selectable (Electric hot wire pilot light for gas jet)  |
| <b>Temperature Probe:</b>                       | Class A Platinum Resistance Thermometer, °C or °F   |
| <b>Flash Detection:</b>                         | Thermocouple  |
| <b>Gas Supply (gas jet ignitor):</b>            | Laboratory gas, 3kPa (0.44psi) maximum pressure   |
| <b>Cooling Supply (Air): (Liquid):</b>          | Dry compressed air, 670kPa (100psi) maximum pressure<br>Water (5°C and above)<br>Water/Ethylene Glycol 50/50 (below 5°C)<br>125kPa (18psi) maximum pressure |
| <b>Voltage:</b>                                 | 110/120V or 220/240V, 50/60Hz, (switchable)   |
| <b>Power:</b>                                   | 1kW   |
| <b>Size (HxWxD):</b>                            | 34 x 47 x 42cm  |
| <b>Weight:</b>                                  | 20kg  |

## ACCESSORIES - Small Scale Closed Cup Module

## OPTIONAL

|           |   |
|-----------|---|
| 00100-0   | SETALERT REMOTE PAGING SYSTEM                 |
| 13770-311 | SYRINGE, 4ml, for tests above 100°C.          |
| 23680-0   | SPIRIT LEVEL                                  |
| 34003-0   | RS232 CABLE, for connecting instrument to PC. |
| 34006-0   | CRYOSTAT (-20°C). For tests to -10°C.         |
| 34008-0   | COOLING BOOST MODULE                          |
| 34010-0   | DRAUGHT SCREEN (recommended for gas tests).   |
| 81002-2   | PRINTER                                       |

## CALIBRATION VERIFICATION &amp; TRAINING

|         |  |
|---------|--|
| 99878-3 | SMALL SCALE CERTIFIED FLASH POINT MATERIAL, (for details see page 156).            |
| 99910-0 | PORTABLE DIGITAL BAROMETER, 750 to 1500 mbar absolute, (for details see page 145). |
| 99928-0 | SETATHERM, (for details see page 129).   |

## CONSUMABLES - Small Scale Closed Cup Module

| 2 Years/2000 Tests |   | Qty |
|--------------------|---|-----|
| 13740-004          | SAMPLE WELL O RING, viton, pack of 5    | 1   |
| 13770-004          | SAMPLE WELL O RING, silicone, pack of 5 | 1   |
| 34200-005          | CO-AXIAL IGNITOR                        | 2   |
| 34700-001          | FLASH DETECTOR                          | 2   |
| 81002-001          | PRINTER PAPER, pack of 20 rolls         | 1   |