

PSL Certified Reference Materials

Certified Reference Materials

PSL Viscosity Reference Standards are suitable for calibration and verification of all types of viscometers. The accuracy of Glass Capillary viscometers, rotational viscometers & falling ball is easily monitored by running a certified viscosity standard. Certification in accordance with UKAS (ISO 17025) guarantees a third party control result - without the inconvenience of having to ship the viscometer for recalibration or verification.

The UKAS (ISO 17025) certification of the PSL reference standards documents direct traceability to the internationally agreed standard value of distilled water at 20 degrees C. At the PSL calibration laboratory a primary viscosity scale is maintained in accordance with internationally agreed procedures.




Each bottle with reference liquid has a unique number and is directly traceable to the primary viscosity scale held in the laboratory. Calibration details on the certified reference materials are obtained using the master viscometer procedures detailed in ASTM D2162.

The accreditation of the PSL Laboratory includes kinematic and dynamic viscosity. Further details available on www.ukas.org, laboratory No. 0247.

All Viscosity Reference Standards from the PSL Calibration Laboratory meet the precision requirements set out in ASTM D445/D446.



PSL Viscosity Reference Standards are supplied in 500 ml sturdy glass bottles complete with calibration certificate and MSDS. The calibration values are also printed on the bottle label. The label is laminated for durability in the laboratory environment.

ACCREDITED TO ISO 17025			
CERTIFICATE OF CALIBRATION			
ISSUED BY THE PSL CALIBRATION LABORATORY			
DATE OF ISSUE	19 December 2008	SERIAL NUMBER	18887
Poulten Selfe & Lee Ltd Russell House, Burnham Business Park, Burnham-on-Crouch, Essex CM0 8TE, England			
		Tel: +44 (0) 1621 787100 Fax: +44 (0) 1621 787175 email: info@rheotek.com www.rheotek.com	
			
Supplied to:	Poulten Selfe and Lee Ltd. Russell House, Burnham-on-Crouch, Essex CM0 8TE		
Labelled:	ASTM reference:	S200	PAGE 1 OF 2 PAGES
	Product code:	2700-V11	APPROVED SIGNATORY
	Batch number:	21518	
	Date of calibration:	11 December 2008	
	Date of expiry:	18 December 2010	
	Bottle number:	55303	 Quality Manager S.P. Sims
A sample taken from this calibrated reference oil was measured against reference viscometers traceable to ASTM D 2162 held at the Laboratory and the following values were determined.			
CERTIFIED VALUES			
Temperature	Kinematic viscosity	Dynamic viscosity	Density
°C	mm ² /s	mPa.s	g/cm ³
20.00	657.5	583.0	0.8867
25.00	458.4	405.1	0.8838
40.00	178.5	156.1	0.8745
50.00	104.7	90.94	0.8686
80.00	30.30	25.77	0.8504
100.00	16.47	13.81	0.8383
UNCERTAINTIES			
Kinematic viscosity	Dynamic viscosity	Uncertainty	
mm ² /s	mPa.s	%	
Up to 6	Up to 6	+/- 0.22	
6 to 19	6 to 19	+/- 0.22	
19 to 100	19 to 100	+/- 0.27	
100 to 450	100 to 450	+/- 0.34	
451 to 1750	450 to 1750	+/- 0.46	
1750 and over	1750 and over	+/- 0.70	
These uncertainties are largely random in nature.			
The reported expanded uncertainty is based on a standard uncertainty multiplied by a coverage factor k=2, providing a level of confidence of approximately 95%. The uncertainty evaluation has been carried out in accordance with UKAS requirements.			

Sample calibration certificate

Catalogue #	Viscosity Standard	Nominal Values	Kinematic Viscosity (KV) in cSt, mm ² /s, Dynamic Viscosity (DV) in cP, mPa.s, Saybolt Viscosity (SUS) values						
			20°C	25°C	37.78°C	40°C	50°C	98.89°C	100.00°C
			68.00 °F	77.00 °F	100.00 °F	104.00 °F	122.00 °F	210.00 °F	212.00 °F
2700-V01H	N.4	KV	0.47	0.45	0.41	0.40			
		DV	0.31	0.29	0.27	0.26			
2700-V02H	N.8	KV	1	0.98	0.76	0.75			
		DV	0.77	0.72	0.57	0.56			
2700-V03H	N1.0	KV	1.3	1.2	0.98	0.97			
		DV	1	0.93	0.77	0.76			
2700-V04	N2	KV	2.9	2.6	2.1	2	1.7		
		DV	2.2	2	1.6	1.5	1.3		
2700-V05	S3	KV	5	4.4	3	2.9	2.6	1.3	1.3
		DV	4.1	3.6	2.5	2.4	2.1	0.98	0.98
2700-V06	S6	KV	11	8.9	5.8	5.7	4.6	1.9	1.9
		DV	8.8	7.4	4.9	4.8	3.7	1.5	1.5
2700-V07	N10	KV	21	17	11	10	7.5	2.7	2.7
		DV	17	14	10	9	6.2	2.2	2.1
2700-V08	S20	KV	47	37	20	18	13	4.4	4
		DV	40	31	17	16	11	3.3	3.2
		SUS			100				
2700-V09	N35	KV	95	72	38	32	23	5.9	5.8
		DV	82	62	33	27	19	4.8	4.7
		SUS			170				
2700-V10	S60	KV	160	120	61	54	35	7.8	7.7
		DV	140	104	54	47	30	6.4	6.3
		SUS			280				
2700-V11	N100	KV	320	230	110	97	59	11	11
		DV	280	200	94	84	51	9.3	9.1
		SUS			500				
2700-V12	N140	KV	400	300	160	140	90	19	18
		DV	350	260	135	120	77	16	15
2700-V13	S200	KV	660	460	210	180	110	18	17
		DV	590	410	180	150	91	15	14
		SUS			930			86	
2700-V14	N230	KV	860	600	270	230	145	23	21
		DV	770	535	230	190	120	19	17
2700-V15	N350	KV	1400	920	370	310	180	25	24
		DV	1200	790	320	270	150	21	20
		SUS						110	
2700-V16	N415	KV	1900	1240	500	415	240	35	34
		DV	1630	1065	430	360	200	29	28
2700-V17	S600	KV	2400	1600	560	520	290	37	35
		DV	2100	1400	490	450	240	31	29
		SUS						150	130
2700-V18	N730	KV	3390	2260	790	730	410	52	49
		DV	2970	1980	690	630	340	43	40
2700-V19	N1000	KV	4800	3100	980	940	520	58	55
		DV	4100	2700	840	800	450	48	45
2700-V20	N1300	KV	6760	4365	1380	1320	730	81	77
		DV	5775	3800	1180	1120	630	67	63
2700-V21	S2000	KV	8600	5600	2000	1700	880	85	81
		DV	7500	4800	1750	1500	760	70	68
		SUS						361	
2700-V22	N4000	KV	18000	11000	3300	3400	1700	140	130
		DV	16000	10000	2800	2900	1500	120	112
2700-V23	S8000	KV	35000	22000	6400	6700	3200	240	220
		DV	31000	20000	6200	5900	2700	210	190
2700-V24	N15000	KV	65000	41000	17000	13000	5800	400	370
		DV	58000	37000	15000	11000	5100	350	320
2700-V25	S30000	KV		82000	32000	23000	11000		670
		DV		74000	30000	21000	9900	640	580