



## King Refrigeration Cloud and Pour Point Baths

### Floor Model Units

Model Number	Number of Baths	Temperatures Degrees C					Sleeves per Bath	Approx. Size (W x L x H)
CP600								Custom design by request
CP601	1						4	Adjustable: -70 through 0 23 x 32 x 36
CP602	1						4	Adjustable: -51 through 0 23 x 20 x 36
CP603	1						4	Adjustable: -34 through 0 23 x 20 x 36
CP604	3	-70	-51	-34			4	23 x 42 x 36
CP605	3		-51	-34	-18		4	23 x 32 x 36
CP606	3			-34	-18	0	4	23 x 30 x 36
CP607	3			-34	-18	0	8	23 x 32 x 36
CP608	4		-51	-34	-18	0	4	23 x 42 x 36
CP609	4		-51	-34	-18	0	8	23 x 44 x 36

### Table Top Model Units

CP610	1						4	Adjustable: -70 through 0 21 x 21 x 22
CP611	1						4	Adjustable: -51 through 0 21 x 21 x 22
CP612	1						4	Adjustable: -34 through 0 23 x 20 x 26
CP613	3		-51	-34	-18		4	23 x 32 x 26
CP614	3			-34	-18	0	4	23 x 30 x 26
CP615	3			-34	-18	0	8	23 x 32 x 26

## Accessories

Part Number	Description
CP620	Test hole Stopper
CP621	Top Insulating Sheet
CP622	Heat Cable
CP623	Felt Sleeve
CP624	Cork, Solid
CP625	Retaining Ring
CP626	Sleeve, Insert
CP627	Foam Insulating Insert
CP628	Test Jar Stopper, Pierced Cork
CP629	Test Jar Stopper, Pierced Cork
CP633	Copper Sleeve
CP634	Cork Ring
CP635	PTC Thermocouple
CP636	Microprocessor (Specify C or F)

### Cloud Point Jars:

**DCP-220-M**

**Cases of 24**

