

## OIL CONTENT OF WAXES

## 20800-0 SETA OIL CONTENT OF WAX APPARATUS

ASTM D721; D3235; IP 158; ISO 2908; NF T60-120;  
FTM 791 5431

- Temperature controlled hot air Evaporation Unit
- Four 15ml bottle capacity
- Rack and pinion operated air tubes
- Three filter stick capacity Cooling Bath
- Kerosine/Solid CO<sub>2</sub> cooling medium

The content of oil and solvents in waxes affect the melting point, hardness and other properties.

The Seta Oil Content of Wax Apparatus comprises a cooling bath for three filter stick assemblies and a temperature controlled hot air evaporation assembly to hold four 15ml weighing bottles.

*Note: Air pump (Seta 14030-2 recommended) required.*



20800-0



SETAPLUS EXTENDED WARRANTY

SPECIFICATIONS	
Voltage:	110/120V, 50/60Hz 220/240V, 50/60Hz
Power:	65W
Heat Lamps:	60W, Tungsten Filament (BC)
Air Flowrate:	2-3 litre/min, per jet
Cooling Bath Capacity:	1.25 litre, total
Cooling Medium:	Kerosine and solid CO <sub>2</sub>
Size (HxWxD):	28 x 27 x 39cm (Cooling Bath) 21 x 21 x 27cm (Evaporation Unit)
Weight:	8.4kg

### CONSUMABLES - Oil Content of Wax Apparatus

2 Years / 2000 Tests		Qty
20480-009	SPRING, pack of 20	1
20850-001	AIR PRESSURE REGULATOR, pack of 2	1
20850-004	WEIGHING BOTTLE, 15ml, pack of 4	1
20850-005	FILTER STICK/COOLING TUBE, pack of 2 pairs	4

### ACCESSORIES - Oil Content of Wax Apparatus

#### REQUIRED

**20850-0** GLASSWARE, comprises air pressure regulator, filter stick and cooling tube with two assembly springs, and four 15ml weighing bottles.

#### THERMOMETERS:

ASTM71F/IP72F; ASTM71C/IP72C

#### OPTIONAL

DRY ICE MAKERS, (for details see page 140).

**14030-2** AIR PUMP, vacuum and pressure, oil-free, (for details see page 136).

**99100-2** ANALYTICAL BALANCE, (for details see page 141).

**45000-0** SETAPLUS EXTENDED 12 MONTH WARRANTY

*Note: Other items which may be required include: Pipettes, Wire stirrer.*